

**WORCESTER POLYTECHNIC INSTITUTE  
HEALTH PHYSICS PROCEDURE HP-19  
GAS-FLOW PROPORTIONAL COUNTER PLATEAU CHECK**

1. PURPOSE:

To ensure that the voltages posted on the gas-flow proportional counter are within the operating plateau.

2. FREQUENCY:

This procedure shall be performed semi-annually.

3. MATERIALS, TOOLS, AND EQUIPMENT:

3.1. Form(s):

3.1.1. Gas Flow Proportional Counter Plateau Check

3.2. Alpha standard button source  
(i.e.: Am-241)

3.3. Beta standard button source  
(i.e.: C-14)

4. PRECAUTIONS:

4.1. Ensure that all health physics practices are followed throughout the survey.

4.2. Utilize the concepts of time, distance, and shielding to maintain exposure as low as reasonably achievable.

5. INSTRUCTIONS:

5.1. Open the large grey valve on the gas tank fully.

5.2. Adjust the small regulating valve so that the gauge reads a gas pressure between 3 and 5 psi.

5.3. Record the gas pressure on the form.

5.4. Ensure that the HV button is in the off position (the button should not be depressed).

- 5.5. Turn on the power, then the HV.
- 5.6. Open the drawer to ensure that the chamber is empty. Close the drawer and turn the knob to the right to check the gas flow. Proper flow rate is approximately 2-3 bubbles/second (the knob may need to be jiggled before the bubbles start to flow).

**NOTE:** Never operate the Gas-Flow Proportional Counter without gas flow.

- 5.7. Determine the alpha plateau
  - 5.7.1. Record the alpha radiation source selected, its reference activity, and its reference date.
  - 5.7.2. Place the source in the chamber.

**NOTE:** Generally, the source is located on the bottom of the button source. Ensure that the source is facing upward.

- 5.7.3. Dial in the first voltage listed on the form.
- 5.7.4. Select a count time of 1 minute (time selector will read 0010).
- 5.7.5. Select a purge time of 144 seconds.
- 5.7.6. Depress reset.
- 5.7.7. Complete a count for each alpha voltage listed on the form (a purge time of 12s may be used for the successive counts).

- 5.8. Repeat step 6 using the beta radiation source and the beta voltages.

**NOTE:** To promote longevity of the machine, dial down the high voltage before depressing the "Add 1000 V" button.

- 5.9. Make two separate graphs, one for the Alpha Data and one for the Beta Data. Plot count rate vs. voltage.
- 5.10. The alpha plateau should appear between 900 - 1200 V and the beta plateau should appear between 1700 - 1800 V. If a relatively flat plateau region does not appear on the graph, take additional counts using smaller voltage increments.
- 5.11. Replace the plateau values listed on the gas-flow proportional counter as necessary. Record the plateau voltages on the form (whether they have changed or not).
- 5.12. Submit the graphs and the form to the RSO for review. Compare the current results to previous results (they should be consistent).
6. RESTORATION:
- None
7. REFERENCES:
- 7.1 Nuclear Measurements Corp. Nuclear Instrumentation Proportional Counting System NMC Model PC-5 Instruction Manual