Worcester Polytechnic Institute AED Inspection & Maintenance Checklist Appendix B

AED Serial Number:	Location:				
	<u>Date</u>	<u>Date</u>	<u>Date</u>	<u>Date</u>	Date
Check ASI flashing green.	<u> </u>	<u> </u>	<u> 2410</u>	<u> </u>	<u> </u>
Check unit and accessories for damage, dirt, and contamination. Clean or replace, as necessary.					
Check battery pack and pads for expiration dates.					
Comments:					
Inspected by (initials):					